

A woman is lying in a hospital bed, wearing a light blue hospital gown. She has a Philips blood pressure cuff on her left arm. The cuff is grey and has the Philips logo and 'Non-Invasive Blood Pressure' printed on it. She is looking towards the camera with a neutral expression. The background shows a hospital room with medical equipment and a white pillow.

PHILIPS

Supplies

Choose both **quality and value**

Match your Efficia devices with Philips supplies for pulse oximetry, ECG, and non-invasive blood pressure measurement

Meeting your needs for **measurement performance** **and affordability**



As a healthcare professional, you need reliable performance from your Philips Efficia devices. Choosing consumables that complement your equipment allows you to get the quality you expect from Philips, at an affordable price. Efficia SpO₂ and ECG, and Value Care non-invasive blood pressure monitoring consumables, are integral parts of Efficia measurement solutions. They work together to provide you clinically reliable, accurate measurements.¹

Efficia sensors, cables, and lead sets, and Value Care cuffs are part of our end-to-end solution for patient care. Designed carefully to meet both your clinical and economic needs, they are worthy of the Philips brand – respected all over the world for quality, reliability, and performance.

¹ These cost-effective sensors and cables have been tested and validated for accuracy on Philips Efficia monitors and defibrillator/monitors (models 863300, 863302, 863304, 866199).



Measure oxygen saturation with ease

Providing continuous, non-invasive measurement of arterial oxygen saturation, our family of durable pulse oximetry sensors are designed for cost-effective reuse. Choose from a variety of different sizes to accommodate children and adults.

Our Efficia reusable SpO₂ sensors are:

- Designed for accuracy
- Easy to clean and disinfect
- Easy to fit and reposition

The products below are all validated for use with Efficia monitors 863300, 863302, and 863304, and Efficia defibrillator/monitors 866199. They are also compatible with any other Philips device that lists them as accessories in the device's Instructions for Use.

	Product number	Description	Patient weight	Application site	Cable length
	989803160611	Efficia Pediatric SpO ₂ Sensor	15 to 50 kg (33 to 110 lbs)	Finger	1.5 m (5 ft)
	989803160621	Efficia Adult SpO ₂ Sensor	> 50 kg (> 110 lbs)	Finger	2 m (6.5 ft)
	989803160631	Efficia Adult SpO ₂ Sensor	> 50 kg (> 110 lbs)	Finger	3 m (10 ft)
	M1196A*	Adult/pediatric clip	> 40 kg (> 88 lbs)	Finger	3 m (10 ft)
	M1196S*	Adult/pediatric clip	> 40 kg (> 88 lbs)	Finger	2 m (6.5 ft)
	M1941A*	SpO ₂ Extension Cable**			2 m (6.5 ft)

* Non-Efficia product which is also recommended for use with Efficia.

**May not be used with the 3 m cabled sensors (M1196A and 989803160631).

Combine ECG solutions to fit the situation at hand

Cover all your basic ECG monitoring needs with our family of Efficia reusable ECG cables and leads. Choose between a one-piece combined trunk cable and lead set, or a two-piece solution (a trunk cable connecting with both 3 or 5 lead sets), for additional value.

Our Efficia ECG solutions feature:

- Reusable, easy-to-clean cables
- Color-coded snap or grabber connectors for easy lead placement
- Shielded ECG wires to reduce artifacts

The products below are all validated for use with Efficia monitors 863300, 863302, and 863304, and Efficia defibrillator/monitors 866199. They are also compatible with any other Philips device that lists them as accessories in the device's Instructions for Use.

	Product number	Description	Color coding	Length
	989803160641	Efficia 3/5-lead ECG patient trunk cable	AAMI/IEC	2.7 m (9 ft)
	989803160651	Efficia 3-lead set, grabber	AAMI	1.0 m (3 ft)
	989803160661	Efficia 3-lead set, grabber	IEC	1.0 m (3 ft)
	989803160671	Efficia 3-lead set, snap	AAMI	1.0 m (3 ft)
	989803160681	Efficia 3-lead set, snap	IEC	1.0 m (3 ft)
	989803160691	Efficia 5-lead set, grabber, limb placement	AAMI	1.0 m/1.6 m (3 ft/5 ft)
	989803160701	Efficia 5-lead set, grabber, limb placement	IEC	1.0 m/1.6 m (3 ft/5 ft)
	989803160711	Efficia 5-lead set, snap, limb placement	AAMI	1.0 m/1.6 m (3 ft/5 ft)

Choose both quality and value



Product number	Description	Color coding	Length
989803160721	Efficia 5-lead set, snap, limb placement	IEC	1.0 m/1.6 m (3 ft/5 ft)
989803160731	Efficia combined one-piece trunk cable/ 3-lead set, grabber	AAMI	2.5 m (8 ft)
989803160741	Efficia combined one-piece trunk cable/ 3-lead set, grabber	IEC	2.5 m (8 ft)
989803160751	Efficia combined one-piece trunk cable/ 3-lead set, snap	AAMI	2.5 m (8 ft)
989803160761	Efficia combined one-piece trunk cable/ 3-lead set, snap	IEC	2.5 m (8 ft)
989803160771	Efficia combined one-piece trunk cable/ 5-lead set, grabber	AAMI	2.5 m (8 ft)
989803160781	Efficia combined one-piece trunk cable/ 5-lead set, grabber	IEC	2.5 m (8 ft)
989803160791	Efficia combined one-piece trunk cable/ 5-lead set, snap	AAMI	2.5 m (8 ft)
989803160801	Efficia combined one-piece trunk cable/ 5-lead set, snap	IEC	2.5 m (8 ft)

Choose cuffs that can stand the pressure

Our Value Care reusable non-invasive blood pressure cuffs feature:

- A robust one-piece design
- Comfortable material

• Color-coded tubing to indicate the size

Our Efficia reusable air hoses feature:

- Metal connectors on both ends to enhance durability
- Easy cleaning, validated with a broad range of disinfectant agents

The products below are all validated for use with Efficia monitors 863300, 863302, and 863304. They are also compatible with any other Philips device that lists them as accessories in the device's Instructions for Use.

	Product number	Description	Limb circumference	Color code
	989803160811	Value Care cuff, Infant	10 to 15 cm (4 to 6 in)	Orange
	989803160821	Value Care cuff, Pediatric	14 to 21.5 cm (5.5 to 8.5 in)	Green
	989803160831	Value Care cuff, Small adult	20.5 to 28 cm (8 to 11 in)	Light blue
	989803160841	Value Care cuff, Adult	27 to 35 cm (10.5 to 14 in)	Navy blue
	989803160851	Value Care cuff, Adult Long	27 to 35 cm (10.5 to 14 in)	Navy blue with edge striping
	989803160861	Value Care cuff, Large adult	34 to 43 cm (13.5 to 17 in)	Burgundy
	989803160881	Efficia adult/pediatric NIBP hose, 1.5 m (5 ft)		
	989803160891	Efficia adult/pediatric NIBP hose, 3 m (10 ft)		

We are continuing to expand our range of Efficia supplies to meet your needs. Please contact your sales representative for more information.

Choose both quality and value

© 2015 Koninklijke Philips N.V. All rights reserved.
Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
(Royal Philips) or their respective owners.

www.philips.com



Printed in The Netherlands.
4522 991 06621 * Jan 2015